



## XFlash® Detector 5030 - 30 mm<sup>2</sup> active area for low beam current applications

Bruker AXS Microanalysis' successful XFlash® Silicon Drift Detector (SDD) series includes a detector with high quantum detection efficiency – the XFlash® Detector 5030. Due to its large active area of 30 mm<sup>2</sup>, this liquid nitrogen-free detector is the instrument of choice for installations on transmission electron microscopes (TEMs) and cold field emission scanning electron microscopes (cold FE-SEMs).

Like other XFlash® SDDs the 5030 has an impressive maximum input count rate capability – 750,000 cps. This is achieved in combination with Bruker's hybrid pulse processing unit. Therefore, the detector delivers extremely fast EDS analysis results.

In comparison to a 10 mm<sup>2</sup> SDD the detector's 30 mm<sup>2</sup> active area results in a threefold increase of count rate detection capability at the same beam current. This makes the XFlash® Detector 5030 ideal for low beam current applications and convenient for the examination of sensitive samples using environmental and variable pressure SEMs as well as TEMs.



Thanks to its optimized electron trap, interference-free analysis is guaranteed even at low excitation energies.

At the same time, the XFlash® 5030 achieves a superb energy resolution of  $\leq 127$  eV (Mn K $\alpha$ ) at 100,000 cps. Nevertheless, the detector's finger is no wider than the finger of the 10 mm<sup>2</sup> XFlash® Detector 5010. Its small dimensions allow reduced distance between detector and sample. The low weight of 2.5 kg minimizes mechanical strain on the TEM or SEM electron column.

The XFlash® Detector 5030 features a passive thermoelectric cooling system, which operates without liquid nitrogen and contains no moving parts. Because of this the detector doesn't affect high resolution microscopes through vibrations or unbalanced mass. The virtually maintenance-free XFlash® 5030 can be used within seconds after being switched on, giving you more time for your applications.

## Technical Data

Energy resolution of 127 eV (Mn K $\alpha$ ) guaranteed at 100,000 cps,  
54 eV C K $\alpha$ , 62 eV F K $\alpha$  (in compliance with ISO 15632 : 2002)  
stable over whole OCR range

Detection from beryllium (4) to americium (95)

Maximum pulse load 750,000 cps

Active area of 30 mm<sup>2</sup>

Optimized electron trap for interference-free analysis in the low energy range


Peltier cooling (no liquid nitrogen or other cooling agents required)

Due to compact design, low weight and vibration-free cooling method no image distortion at the SEM

Compatible with all microscopes

### Comparison XFlash<sup>®</sup> 5030 / Si(Li) Detector

Typical energy resolution at Mn K $\alpha$  vs. input count rate



## ● Bruker AXS Microanalysis

Schwarzschildstr. 12  
12489 Berlin  
Germany  
Tel. +49 (30) 670990-0  
Fax +49 (30) 670990-30  
info-ma@bruker-axs.de  
www.bruker-axs-ma.com

Bruker AXS Microanalysis in:

**Brazil**  
Tel. +55 (11) 5052 5030  
info@bruker.com.br  
www.bruker-axs-ma.com

**Mexico**  
Tel. +52 (55) 5601 2599  
info-axs@bruker.com.mx  
www.bruker-axs-ma.com

**P.R. China**  
Tel. +86 (10) 68486946  
info@bruker-axs.cn  
www.bruker-axs.cn

**Nordic Countries**  
Tel. +46 (8) 54480820  
info@bruker-axs.se  
www.bruker-axs-ma.com

**Japan**  
Tel. +81 (45) 4531960  
info@bruker-axs.jp  
www.bruker-axs.jp

**South Africa**  
Tel. +27 (11) 463 6040  
info@bruker.co.za  
www.bruker.com/za

**Korea**  
Tel. +82 (2) 3476 8600  
info@bruker-axs.co.kr  
www.bruker-axs.co.kr

**USA**  
Tel. +1 (609) 771 4400  
info-ma@bruker-axs.com  
www.bruker-axs-ma.com

Or find your local partner at: [www.bruker-axs-microanalysis.com](http://www.bruker-axs-microanalysis.com)